



USER INSTRUCTIONS

DONNING INSTRUCTIONS: Follow these instructions each time the respirator is worn. Inspect prior to use. Examine all the respirator parts for signs of damage including the two headbands and the nosepiece. The respirator should be disposed of immediately upon observation of damaged or missing parts.



STEP 1

Separate the edges of the respirator to fully open it. Slightly bend the nosepiece to form a curve.



STEP 2

Using your index fingers and thumbs, separate both headbands.




STEP 3

Continue holding the headbands and cup the respirator under your chin with the nosepiece facing upwards.



STEP 4

Pull the headbands over your head and let them rest at the base of your neck, below your ears.



STEP 5

Position the top headband to rest at the crown of your head, above your ears.



STEP 6

Using both hands, firmly press down to form the nosepiece to the bridge of your nose and below the eyes. Adjust to secure all edges and perform a seal check.

REMOVAL INSTRUCTIONS: Pull the bottom headband over your head. Still holding the respirator in position, pull the top headband over your head and remove the respirator.

✓ PERFORMING A SEAL CHECK:

It is necessary to perform a seal check prior to each time of use. To check the respirator seal, cup the respirator with both hands, breathe in forcefully. You should feel pressure inside the respirator. If you feel any air leakage on an inhale or exhale, adjust the nosepiece and headbands until you feel a secure fit with no air leakage. Repeat the seal check again to confirm proper fit. **DO NOT enter a contaminated area if a proper seal cannot be obtained.**

CONTRAINDICATIONS: This device is not designed, sold or intended for use except as indicated. Not for use with beards, or other facial hair that prevents direct contact between the face and sealing surface of the respirator.

TIME USE LIMITATIONS: Respirator may be used until damaged, breathing becomes difficult, or contaminated. Discard after every use. Follow national, state, local, and facility infection control guidance and policies.

⚠ WARNING:

When worn properly, this Surgical N95 Respirator helps reduce exposure to certain airborne particles, but does not eliminate exposure to or the risk of contracting any disease or infection. Before use, the wearer must read and understand the instructions provided. Misuse may result in sickness or death.

REMINDER: Before occupational use of this respirator, a written respiratory protection program must be implemented meeting all the local government requirements. In the United States, employers must comply with OSHA 29 CFR 1910.134 which includes medical evaluation, training, and fit testing.

MATERIALS:

- Outer Splash Layer: Spunbond Polypropylene
- Second Splash Layer: Spunbond Polypropylene
- Filter Layer: Meltblown Polypropylene
- Skin-Touch Layer: Spunbond Polypropylene

DESCRIPTION: The Surgical Respirator is NIOSH approved as an N95 particulate filter respirator. It provides 95% filtration efficiency of 0.3 micron particles. It meets the CDC Guidelines for TB exposure and control.

INTENDED USE: This Surgical N95 Respirator is a single-use disposable device intended to be worn by operating room personnel and other healthcare workers to protect the professional and the patient by reducing the exposure to airborne contaminants, and the transfer of microorganisms, body fluids and particulates.

S - SPECIAL OR CRITICAL USER INSTRUCTIONS

This respirator has been approved as a NIOSH N95 filtering facepiece respirator, for use in healthcare settings, as a Surgical N95 Respirator conforming to recognized standards for biocompatibility, flammability, and fluid resistance.

- Fluid Resistant at 80, 120, and 160 mmHg per ASTM F1862
- Flammability rated Class I per 16 CFR 1610
- Inhalation < 35 mm H₂O
- Exhalation < 25 mm H₂O



Manufactured by ivWatch, LLC at our
FDA registered facility
700 Tech Center Parkway, Suite 300
Newport News, VA 23606 USA

MADE IN USA

To order or learn more, visit: blox.ivWatch.com
If you have any questions, contact: support@ivWatch.com

**Part # BLX-9501
QL-1002006 Rev.06**